

# CASSINI<sup>TM</sup>

C-900M

MICROSLICER USE

INSERT SPECIMENS YOU WISH TO STUDY INTO ONE OF THE HOLES OF THE MICRO-SLICER . ROTATE THE KNOB TO CUT YOUR SPECIMEN INTO THIN SLICES. THE MICROSLICER IS AN IDEAL TOOL IN THE MAKING OF SECTION SLIDES.

**CAUTION: THE BLADE OF THE MICROSLICER IS VERY SHARP. HANDLE THE MICRO-SLICER WITH CARE.**



**1**



**2**



**3**